




From Chaos to Control:

A Salivary Diagnostics-Driven Approach to Managing High-Risk Caries Patients



Speaker

Dr. Kim Kutsch



Speaker

Dr. Doug Thompson



Speaker


Dr. Doug Young

In this webinar, you'll learn:

- How salivary diagnostics guide treatment decisions
- Strategies for targeted biofilm management
- Practical steps to stabilize high-caries-risk patients long-term

Stay tuned for a **special offer**
at the **end** of tonight's webinar!

Machiulskiene V, Campus g, Carvalho JC, Dige J, Ekstrand KR et al. Terminology of dental caries and dental caries management: consensus report by ORCA and IADR. Caries Research 2020. 54(1):7-14.





“Dental caries is a biofilm-mediated, diet-modulated, multifactorial, non-communicable dynamic disease resulting in net mineral loss of dental hard tissues.”

Bradshaw DJ, McKee AS,
Marsh PD. Effects of
carbohydrate pulses and pH
on population shifts within
oral microbial communities in
vitro. J Dent Res. 1989
Sep;68(9):1298-302.




“Analysis of the data strongly
suggests that the pH
generated from carbohydrate
metabolism, rather than
carbohydrate availability per
se, is responsible for the
widely reported shifts in
composition and metabolism
of the oral microflora in vivo.”


WHY IS pH SO IMPORTANT?




STEPHAN CURVE (Normal)




STEPHAN CURVE (Snacking)



STEPHAN CURVE (Hyposalivation)



STEPHAN CURVE (S+H)



CRA Form

CRA Form Adults and Children 6+ Four-Step Guide to Cavity Prevention

Name _____

Date _____


1 Risk & Disease Discovery Complete the section below by circling the most appropriate answer for each question.


RISK FACTORS		
SALIVA		
Do you take medications daily? If so, how many?	NO	YES (circled)
Do you feel as though you have a dry mouth at any time of the day or night?	NO	YES
DIET		
Do you drink liquids other than water more than 2 times daily between meals?	NO	YES
Do you snack daily between meals?	NO	YES
BIOFILM		
Do you notice plaque build-up on your teeth between brushings?	NO	YES
DISEASE INDICATORS		
New/Progressing visible cavitations?	NO	YES (circled)
New/Progressing approximal radiographic radioluscencies?	NO	YES (circled)
New/Active white spot lesions?	NO	YES (circled)
Is decay history a concern?	NO	YES (circled)

2 Risk Identification Determine risk based on answers above; one blue response indicates moderate risk, one red response indicates high risk.


Healthy	+ Risk Factors	+ Disease Indicators
1 - Low Risk	2 - Moderate Risk	3 - High Risk (circled)
CDT Code D0601	CDT Code D0602	CDT Code D0603

Once risk level diagnosis has been made, refer to next page for treatment options and reassessment protocol.






Modified from Kornman KS: Clin Infect Dis 28: 520, 1999



Modified from Kornman KS: Clin Infect Dis 28: 520, 1999

Etiology: The Microbiome of Caries

- **Highly acidic** and **aciduric** species associated with caries include *Streptococcus mutans*, *Lactobacillus*, *Actinomyces*, *Bifidobacterium*, and *Scardovia* species.
- Excess acid production leads to a **dysbiotic shift** in the biofilm composition → demineralization of tooth structure → cavities.



Scanning electron micrograph showing co-adhesion of *Streptococcus mutans* with *Candida albicans*
Image acquired by Wanessa Fernandes Matias Regis

- The role of bacteria in the etiology of dental caries is long established, while the **role of fungi has only recently gained more attention.**
- There was a **higher ratio of fungi to bacteria in dentin-involved lesions compared to less severe lesions**
 - *Candida albicans*, *C. dubliniensis*, and *C. tropicalis*.
- *Candida albicans* and *Streptococcus mutans* frequently coexist within dental biofilms and **enhance each other's pathogenic potential.**
- These mixed-species biofilms exhibit **increased acid production** and **greater resistance to antimicrobial agents** compared to single-species biofilms.

Matias Regis, W. F., et al. (2025). Insights into the role of *Streptococcus mutans* and *Candida albicans* in dental biofilm formation and cariogenicity: A literature review. *Cureus*, 17(6), e86159. <https://doi.org/10.7759/cureus.86159>

Spatafora, G., et al (2024). The evolving microbiome of dental caries. *Microorganisms*, 12(1), 121. <https://doi.org/10.3390/microorganisms12010121>

“Are we treating teeth... or are we treating disease?”

Research supports the disease of caries is biofilm mediated and host modulated and in addition to causing carious lesions in the teeth it affects:

- Serum hsCRP
- Stroke risk (up to three-fold)
- Advancement of cardiovascular disease and atherosclerotic plaque
- Microbleeds in the brain
- Nitric Oxide production impairment
- Inflammatory bowel disease
- Tumors
- Infective endocarditis
- Etc.





FIGURE 1
Multi-organ damage caused by *S. mutans* in systemic diseases.

AlShammari, A., et al. (2024). The association between dental caries and serum CRP in the US adult population: Evidence from NHANES 2015–2018. *BMC Public Health*, 24, 2210. <https://doi.org/10.1186/s12889-024-19681-6>

Fang, Y., et al. (2024). Roles of *Streptococcus mutans* in human health: Beyond dental caries. *Frontiers in Microbiology*, 15, 1503657. <https://doi.org/10.3389/fmicb.2024.1503657>

Management Models

- **Traditional management** is based on the **surgical model: cutting and restoring the tooth to remove caries lesions.**
 - This restorative cycle can be **short-sighted** by not addressing full etiological source of caries.
- **The medical model of caries management is more conservative, preventive, and patient-centered.**
- The medical model involves four main principles
 1. **Control of bacterial infection**
 2. **Reduction of risk levels**
 3. **Remineralization of teeth**
 4. **Long-term follow-up.**



Caries Management | Time for A New Management Model? | Control the Bacterial Infection

Sample Report
 Date of Birth: 08/17/2000 (22 yrs)
 Gender: Female
 Patient ID: 920-C
 Patient Location: Test Site A

Ordering Provider
 Ronald McGlennen MD
 PACO Primary Carel Drive
 Suite 150
 Eden Prairie, MN 55544
 855-672-5362

Specimen#: 560020029
 Accession#: 20220-CR38
 Specimen: Oral Rinse/P1
 Collected: 06/11/2023
 Received: 06/12/2023 15:14
 Reported: 06/13/2023 10:24

MOLECULAR DETECTION OF CARIES RISK PATHOGENS

Reason for Testing: Screening/Risk Assessment
 Dental History: Recent/Current History of Caries
 Related Info: Not Provided

CARIES RISK PATHOGENS

Bacteria	Test Result
Streptococcus mutans	Detected
Streptococcus sobrinus	Not Detected
Lactocaseibacillus casei	Not Detected

SUMMARY OF RESULTS

Interpretation:
 S. mutans has been detected in this sample. The combined amounts of these pathogens indicates a HIGH RISK of the development and progression of dental caries.

Significance:
 The detection of S. mutans DNA indicates presence of these organisms in the oral cavity. This is an indicator of a change in the oral microbiome, which typically leads to the creation of an acidic environment, a contributing factor in caries formation.

Methodology: Genomic DNA is extracted from the submitted sample and tested for 3 species of bacteria (Streptococcus mutans, Streptococcus sobrinus, Lactocaseibacillus casei). Bacterial DNA is assayed using qPCR (TM) based oligonucleotides and real-time quantitative polymerase chain reaction (qPCR). Bacterial levels are represented in genome copies per mL of sample. Risk ranges were derived from patient testing: Low risk (0 copies/mL), Moderate risk (only to 10,000 copies/mL), High risk (greater than 10,000 copies/mL). The analytical and performance characteristics of this laboratory-developed test (LDT) was determined by OralDNA.

OralDNA OraRisk Caries
 Tests for *S. mutans*, *S. sobrinus*, and *L. casei*
 Provides logarithmic amount, if detected

FIDALAB Comprehensive Cariogenic Pathogens Test LAB RESULTS

Patient Information: Last Name: Romzick, First Name: Peter, Date of Birth: 04/24/1959, Gender: Male
 Ordering Provider: Integrative Oral Medicine
 Sample Information: Accession: F36405-CCAR, Specimen type: Oral Fluid, Collected Date: 3/17/2025, Receiving Date: 3/20/2025

A. Bacteria Pathogens Results

Type	Organism	Copyes/ml (range)	10 ⁷	10 ⁶	10 ⁵	10 ⁴	10 ³	10 ²	10 ¹
A	Streptococcus mutans	52,413	█						
	Streptococcus sobrinus	N/A							
	Lactocaseibacillus casei	737,284	█	█	█				
B	Lactobacillus fermentor	1435,154	█	█	█	█			
	Lactobacillus rhamnosus	1335,157	█	█	█	█			
	Bifidobacterium dentium	1038,053	█	█	█	█			
	Streptococcus oralis	N/A							
C	Actinomyces viscosus	4929,965	█	█	█	█	█		
	Weillonia spp.	37408,757	█	█	█	█	█	█	
	Saurelavisceae spp.	24,009	█	█	█	█			
	Propionibacterium acidiflavum	N/A							

N/A Below Detection Limit * << 10³ Copies/sample Bacterial pathogens composite score:301.12

B. Fungal Pathogens Results

Candida albicans	Detected
Candida dubliniensis	Not Detected
Candida parapsilosis	Not Detected
Candida glabrata	Not Detected
Candida krusei	Not Detected
Candida tropicalis	Not Detected
Candida kefyr	Not Detected
Candida guilliermondii	Not Detected
Candida lusitanae	Not Detected
Candida rugosa	Not Detected

FidaLab Comprehensive Caries
 Tests for *S. mutans*, *S. sobrinus*, and nine more bacteria
 Tests for 10 *Candida* spp.
 Provides logarithmic amount and composite score

SIMPLY PERIO™ SPECTRUM SIMPLY TEST™

FINAL REPORT
 PERIODONTAL PREMIUM
 Sample Type: Saliva
 Reported: 11/21/2023 17:14

PATIENT INFO: First18 Last18, DOB: 01/18/2000, SEX: MALE
 SAMPLE INFO: Specimen#: ST01810, Collected: 11-21-2023 14:20, Received: 11/21/2023 15:19
 ORDERING PROVIDER: John Doe MD, NPI: 1234567890, PHONE: [REDACTED]

LOW MEDIUM HIGH NON-VIRAL TARGETS (In Range Results)

All displayed values are in genomic copies x1000/ml, except Fusobacterium nucleatum which is in genomic copies x10,000/ml.

Spectrum Diagnostics SimplyPerio
 Tests for *S. mutans*, *S. sanguinis*, and *C. albicans*
 Provides logarithmic amount, if detected



CariScreen Meter
 “CariScreen score had a strong positive correlation $r = 0.76$ with total cell count, a positive correlation with *MS* counts $r = 0.69$, and a positive correlation to caries risk status $r = 0.55$ with high significance $p = 0.000001$.”

Case Study: Monica L.
Dysbiotic Microbiome

Dental History Review

- In 2000 at age 20, had **20 cavities** which were taken care of at University of Michigan Dental School along with orthodontia
- **Patient stated she still gets constant cavities "...no matter much I brush or floss or Waterpik..."**
- In 2010 at age 30, had **all amalgam fillings removed** and replaced.
- Experiences **unpleasant tastes and odors** from her mouth
- due to **food getting stuck**
- Reports **extreme gum recession** and that she has been told she has **very short roots** from other dentists
- Reports **constant cavities no matter how good her oral hygiene habits are or how frequently she gets her teeth cleaned**
- **Sought a second opinion after being advised to extract all teeth**

How would you rate the condition of your mouth?

Excellent	No	Fair	No
Good	No	Poor	YES

Previous Dentist?
Shoha

How long have you been a patient? (months/years)
2017

Date of most recent dental exam?
within 6 months

Date of most recent x-rays?
within 6 months

Date of most recent treatment (other than a cleaning)
within 6 months

I routinely see my dentist every

3 mo	No	6 mo	YES	Not routinely	No
4 mo	No	12 mo	No		

What is your immediate concern?
Decay

PLEASE ANSWER YES OR NO TO THE FOLLOWING:

PERSONAL HISTORY

1. Are you fearful of dental treatment? How fearful, on a scale of 1 (least) to 10 (most)	YES	5
2. Have you had an unfavorable dental experience?	YES	plenty
3. Have you ever had complications from past dental treatment?	YES	didn't heal from root canal
4. Have you ever had trouble getting numb or had any reactions to local anesthetic?	No	
5. Did you ever have braces, orthodontic treatment or had your bite adjusted, and at what age?	YES	Lots of orthodontia as teen
6. Have you had any teeth removed, missing teeth that never developed or lost teeth due to injury or facial trauma?	YES	Removed tooth in right rear upper due to not wanting to give it a root canal

BITE AND JAW JOINT








21. Do you have problems with your jaw joint? (pain, sounds, limited opening, locking, popping)	No
22. Do you feel like your lower jaw is being pushed back when you try to bite your back teeth together?	No
23. Do you avoid or have difficulty chewing gum, carrots, nuts, bagels, baguettes, protein bars, or other hard, dry foods?	No
24. In the past 5 years, have your teeth changed (become shorter, thinner, or worn) or has your bite changed?	No
25. Are your teeth becoming more crooked, crowded, or overlapped?	No
26. Are your teeth developing spaces or becoming more loose?	No
27. Do you have trouble finding your bite, or need to squeeze, tap your teeth together, or shift your jaw to make your teeth fit together?	No
28. Do you place your tongue between your teeth or close your teeth against your tongue?	No
29. Do you chew ice, bite your nails, use your teeth to hold objects, or have any other oral habits?	No
30. Do you clench or grind your teeth together in the daytime or make them sore?	No
31. Do you have any problems with sleep (i.e. restlessness or teeth grinding), wake up with a headache or an awareness of your teeth?	No
32. Do you wear or have you ever worn a bite appliance?	No

9. Have you ever noticed an unpleasant taste or odor in your mouth?	YES	I realize now that it's because I had food stuck where I didn't know about it
10. Is there anyone with a history of periodontal disease in your family?	No	
11. Have you ever experienced gum recession, or can you see more of the roots of your teeth?	YES	Extreme
12. Have you ever had any teeth become loose on their own (without an injury), or do you have difficulty eating an apple?	YES	Shoha said only one tooth was loose but other dentists have said very short roots
13. Have you experienced a burning or painful sensation in your mouth not related to your teeth?	No	

TOOTH STRUCTURE


14. Have you had any cavities within the past 3 years?	YES	Constant cavities no matter how good my oral hygiene is or if I get cleanings ever 3 months instead of 6 I feel I have a lot of saliva
15. Does the amount of saliva in your mouth seem too little or do you have difficulty swallowing any food?	YES	
16. Do you feel or notice any holes (i.e. pitting, craters) on the biting surface of your teeth?	No	
17. Are any teeth sensitive to hot, cold, biting, sweets, or do you avoid brushing any part of your mouth?	No	
18. Do you have grooves or notches on your teeth near the gum line?	No	
19. Have you ever broken teeth, chipped teeth, or had a toothache or cracked filling?	No	
20. Do you frequently get food caught between any teeth?	YES	In gum recession

FMX | 11-06-2023



November 2023

- Consult for OralDNA MyPerioPath + Candida and FidaLab salivary diagnostic tests.



MOLECULAR IDENTIFICATION OF CANDIDA SPECIES IN THE OROPHARYNX

Test Results
Candida Species: Negative

Signs and Symptoms of Oral Candidiasis

- Often no symptoms
- "Burning Mouth Syndrome"
- Metallic or acidic or salty taste

Causes

- Various Candida species, most often C.albicans
- Underlying systemic disease
- Immunosuppression

Interpretation:
This sample is negative for DNA of the following Candida species: C. albicans, C. glabrata, C. krusei, C. parapsilosis, C. tropicalis, C. kefyr, C. guilliermondii, C. lusitanae, and C. rugosa. See comments.

Comments:

Significance: Oral candidiasis (Oral thrush) is caused by a candida fungus that can infect the mouth. Oral thrush can cause creamy white lesions, usually on the tongue or inner cheeks, and may spread to the roof of mouth, gums, tonsils, or the back of the throat. It is not a type of infection that can be passed on to others.

Risk: Oral thrush most commonly affects people who wear dentures. People who have difficulties keeping their mouth clean, people with diabetes and those who take steroids are also at a higher risk of developing the condition. Some antibiotics may cause thrush. Certain antibiotics encourage the infection to recur, especially if taken over a long period of time. Very rarely, oral thrush may be one of the early signs of HIV.

		Results						
Type	Organism	Copies/oral rinse sample (x10 ⁶)	10 ³	10 ⁵	10 ⁷	10 ⁸	10 ⁹	
A	Streptococcus mutans	314.426	[Red bar]					
	Streptococcus sobrinus	N/A	[None]					
B	Lactobacillus casei	05.793	[Red bar]					
	Lactobacillus fermentum	04.510	[Red bar]					
	Lactobacillus rhamnosus	00.301	[Red bar]					
	Bifidobacterium dentium	287.053	[Red bar]					
C	Stactia exigua	N/A	[None]					
	Actinomyces viscosus	4840.273	[Yellow bar]					
	Veillonella spp	8629.421	[Yellow bar]					
	Scardovia wiggsiae	05.600	[Yellow bar]					
	Propionibacterium acidifaciens	64.059	[Yellow bar]					

N/A Below Detection Limit * < 10³ Copies/sample **Cariogenic Pathogen Composite Score: 231.85**

Interpretation

- Streptococcus mutans was detected.
- Streptococcus sobrinus was not detected.
- The patient has **High risk (Composite Score: 231.85)** for dental caries (tooth decay).

Comments

- Over 90% of U.S. adults have dental caries in permanent teeth, and about 13% of children and adolescents have untreated dental caries in permanent teeth.
- The first step in dental caries (also known as tooth decay) development is the overgrowth of cariogenic bacteria in the oral cavity, which in turn generates acid that breaks down the tooth enamel.
- Proper oral health including professional cleaning to remove dental plaques can reduce the risk for dental caries development.
- Avoid foods and beverages with high sugar content that are acidic (most beverages other than water are acidic).
- Prescription strength fluoride used daily can reduce the risk for dental caries.

OralDNA MyPerioPath

- No high-risk red-complex pathogens
- Fn only present orange-complex pathogen
- Cs present

OralDNA OraRisk Candida

- Negative for *Candida* spp.
C. albicans, *C. glabrata*, *C. krusei*, *C. parapsilosis*, *C. tropicalis*, *C. kefyr*, *C. guilliermondii*, *C. lusitanae*, and *C. rugosa*

FidaLab Cariogenic Pathogen Test

- Composite score of 231.85 - **HIGH CARIES RISK**
- S. mutans* present - **314.426 x 10⁷**

November 2023 – February 2024

- **CHIEF CONCERN: CARIES**
- Recommended **CariFree Treatment Rinse and gel toothpaste.**
- **Hopeless teeth removed and Zirconia implants scheduled for early June.**
- Continue **CariFree Treatment Rinse** until profile is satisfactory.
- **Restorative work on non-hopeless teeth completed** at Integrative Oral Medicine.
- FidaLab Caries **retest scheduled for every three to six months.**



November 2023 - December 2024

Extractions, implants, restorations to restore function. CariFree Treatment Rinse and 3 month recare appointments the entire time.

November 2023

Results								
Type	Organism	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	
A	Streptococcus mutans	314.426	[Red bar]					
	Streptococcus sobrinus	N/A						
B	Lactobacillus casei	08.793	[Red bar]					
	Lactobacillus fermentum	04.510	[Red bar]					
	Lactobacillus rhamnosus	00.301	*					
	Bifidobacterium dentium	287.053	[Red bar]					
	Slackia exigua	N/A						
C	Actinomyces viscosus	4840.273	[Yellow bar]					
	Veillonella spp	8629.421	[Yellow bar]					
	Scardovia wiggisiae	08.600	[Yellow bar]					
	Propionibacterium acidifaciens	64.059	[Yellow bar]					

N/A Below Detection Limit * <- 10³ Copies/sample **Cariogenic Pathogen Composite Score:231.85**

Interpretation

- Streptococcus mutans was detected.
- Streptococcus sobrinus was not detected.
- The patient has **High risk (Composite Score: 231.85)** for dental caries (tooth decay).

Comments

- Over 90% of U.S. adults have dental caries in permanent teeth, and about 13% of children and adolescents have untreated dental caries in permanent teeth.
- The first step in dental caries (also known as tooth decay) development is the overgrowth of cariogenic bacteria in the oral cavity, which in turn generates acid that breaks down the tooth enamel.
- Proper oral health including professional cleaning to remove dental plaques can reduce the risk for dental caries development.
- Avoid foods and beverages with high sugar content that are acidic (most beverages other than water are acidic).
- Prescription strength fluoride used daily can reduce the risk for dental caries.

March 2024

Results								
Type	Organism	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	
A	Streptococcus mutans	141.956	[Red bar]					
	Streptococcus sobrinus	N/A						
B	Lactobacillus casei	02.994	[Red bar]					
	Lactobacillus fermentum	11.984	[Red bar]					
	Lactobacillus rhamnosus	24.338	[Red bar]					
	Bifidobacterium dentium	93.526	[Red bar]					
	Slackia exigua	N/A						
C	Actinomyces viscosus	397.133	[Yellow bar]					
	Veillonella spp	11222.268	[Yellow bar]					
	Scardovia wiggisiae	03.005	[Yellow bar]					
	Propionibacterium acidifaciens	22.826	[Yellow bar]					

N/A Below Detection Limit * <- 10³ Copies/sample **Cariogenic Pathogen Composite Score:106.87**

Interpretation

- Streptococcus mutans was detected.
- Streptococcus sobrinus was not detected.
- The patient has **High risk (Composite Score: 106.87)** for dental caries (tooth decay).

Comments

- Over 90% of U.S. adults have dental caries in permanent teeth, and about 13% of children and adolescents have untreated dental caries in permanent teeth.
- The first step in dental caries (also known as tooth decay) development is the overgrowth of cariogenic bacteria in the oral cavity, which in turn generates acid that breaks down the tooth enamel.
- Proper oral health including professional cleaning to remove dental plaques can reduce the risk for dental caries development.
- Avoid foods and beverages with high sugar content that are acidic (most beverages other than water are acidic).
- Prescription strength fluoride used daily can reduce the risk for dental caries.

December 2024

Results								
Type	Organism	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	
A	Streptococcus mutans	00.033	*					
	Streptococcus sobrinus	N/A						
B	Lactobacillus casei	316.729	[Red bar]					
	Lactobacillus fermentum	79.356	[Red bar]					
	Lactobacillus rhamnosus	76.175	[Red bar]					
	Bifidobacterium dentium	227.077	[Red bar]					
	Slackia exigua	N/A						
C	Actinomyces viscosus	397.176	[Yellow bar]					
	Veillonella spp	1197.225	[Yellow bar]					
	Scardovia wiggisiae	N/A						
	Propionibacterium acidifaciens	00.021	*					

N/A Below Detection Limit * <- 10³ Copies/sample **Cariogenic Pathogen Composite Score:36.03**

Interpretation

- Streptococcus mutans was detected.
- Streptococcus sobrinus was not detected.
- The patient has **High risk (Composite Score: 36.03)** for dental caries (tooth decay).

Comments

- Over 90% of U.S. adults have dental caries in permanent teeth, and about 13% of children and adolescents have untreated dental caries in permanent teeth.
- The first step in dental caries (also known as tooth decay) development is the overgrowth of cariogenic bacteria in the oral cavity, which in turn generates acid that breaks down the tooth enamel.
- Proper oral health including professional cleaning to remove dental plaques can reduce the risk for dental caries development.
- Avoid foods and beverages with high sugar content that are acidic (most beverages other than water are acidic).
- Prescription strength fluoride used daily can reduce the risk for dental caries.

FidaLab Cariogenic Pathogen Test

- Composite score of 231.85 - **HIGH CARIES RISK**
- S. mutans* present - **314.426 x 10⁷**

FidaLab Cariogenic Pathogen Test

- Composite score of 231.85 > 106.87 - **Pt is still considered a HIGH CARIES RISK**
- S. mutans* present - **314.426 x 10⁷ > 141.956 x 10⁷**

FidaLab Cariogenic Pathogen Test

- Composite score of 106.87 > 36.03 - **Pt is still considered a HIGH CARIES RISK but composite score has DECREASED by 84%**
- S. mutans* NOT PRESENT**

January 2025

- **Implant support zirconia crowns on #3, 12, 18, 19, 30, 31.**
- **Agreed to 6-month radiograph recare to monitor for new lesions. Significantly reduced caries risk and new biofilm profile.**
- New FidaLab Cariogenic Pathogen Test and Fidalab Comprehensive Periodontal Pathogen Test administered.



January 13, 2025

- Caries composite score reduction over time: **231.85 > 106.87 > 36.03 > 00.92**

A. Bacterial Pathogens Results										
Type	Organism	Copies/oral rinse sample (x10 ⁶)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰	10 ¹¹	
A	Porphyromonas gingivalis(PO)	N/A								
	Tannerella forsythia(TF)	N/A								
	Treponema denticola(TD)	N/A								
	Prevotella nigrescens(PN)	03.205								
B	Porphyromonas endodontalis(PE)	N/A								
	Prevotella intermedia(PI)	N/A								
	Eubacterium nodosum(EN)	N/A								
	Factor albicans(FA)	N/A								
	Fretobacterium sp. HOT380(FR)	N/A								
	Treponema soerenga(TS)	N/A								
	C	Bacteroides [O-2] sp. HOT274(Ba)	N/A							
		Campylobacter rectus(CR)	00.047	*						
		Parvimonas micra(PM)	00.082	*						
		Desulfobulbus sp. HOT041(De)	N/A							
Eikenella corrodens(EC)		110.025								
Selenomonas sputigena(SS)		110.575								
Eubacterium saeponum(ES)		N/A								
TM7 [G-1] sp. HOT349(TM7)		410.841								
Fretobacterium fastidiosum(FR)		N/A								
Treponema leathemolyticum(TL)		N/A								
D	Treponema maingilium(TM)	N/A								
	Fusobacterium nucleatum ss animalis(FNa)	00.703	*							
	Fusobacterium nucleatum ss polymorphum(FNp)	58.553								
	Treponema sp. HT23(TT)	N/A								
	Fusobacterium nucleatum ss vincentii(FNv)	01.454								
E	Capnocytophaga gingivalis(CG)	37.528								
	Capnocytophaga ochracea(CO)	24.147								
	Capnocytophaga sputigena(CS)	01.728								
E	Aggregatibacter actinomycetemcomitans(AA)	N/A								

FidaLab Perio - Bacteria

- No red-complex bacteria
- Composite score: **00.77**

B. Fungal Pathogens Results	
<i>Candida albicans</i>	Not Detected
<i>Candida dubliniensis</i>	Not Detected
<i>Candida parapsilosis</i>	Not Detected
<i>Candida glabrata</i>	Not Detected
<i>Candida krusei</i>	Not Detected
<i>Candida tropicalis</i>	Not Detected
<i>Candida kefir</i>	Not Detected
<i>Candida guilliermondii</i>	Not Detected
<i>Candida lusitanae</i>	Not Detected
<i>Candida rugosa</i>	Not Detected

C. Viral Pathogens Results	
HHV-1 (HSV1)	Not Detected
HHV-2 (HSV2)	Not Detected
HHV-3 (VZV)	Not Detected
HHV-4 (EBV)	Not Detected
HHV-5 (CMV)	Not Detected
HHV-6	Not Detected
HHV-8	Not Detected

FidaLab Perio – Yeast / Viruses

- **Negative for 10 *Candida* spp.**
- **Negative for HHV-6**

Results							
Type	Organism	Copies/oral rinse sample (x10 ⁶)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹
A	Streptococcus mutans	00.346	*				
	Streptococcus sobrinus	N/A					
B	Lactobacillus casei	02.205					
	Lactobacillus fermentum	06.270					
	Lactobacillus rhamnosus	00.118	*				
	Bifidobacterium dentium	00.027	*				
C	Slackia exigua	N/A					
	Actinomyces viscosus	59.422					
	Veillonella spp	5458.099					
	Scardovia wiggsiae	N/A					
	Propionibacterium acidifaciens	N/A					

N/A Below Detection Limit * <- 10⁵ Copies/sample Cariogenic Pathogen Composite Score:00.92


Interpretation	
• Streptococcus mutans was detected.	
• Streptococcus sobrinus was not detected.	
• The patient has Low risk (Composite Score: 00.92) for dental caries (tooth decay).	

Comments	
• Over 90% of U.S. adults have dental caries in permanent teeth, and about 13% of children and adolescents have untreated dental caries in permanent teeth.	
• The first step in dental caries (also known as tooth decay) development is the overgrowth of cariogenic bacteria in the oral cavity, which in turn generates acid that breaks down the tooth enamel.	
• Proper oral health including professional cleaning to remove dental plaques can reduce the risk for dental caries development.	
• Avoid foods and beverages with high sugar content that are acidic (most beverages other than water are acidic).	
• Prescription strength fluoride used daily can reduce the risk for dental caries.	

FidaLab Cariogenic Pathogen Test


- Composite score of **231.85 > 106.87 > 36.03 > 00.92**
- **Pt is no longer considered HIGH CARIES RISK composite score has reduced significantly**
- ***S. mutans* NOT PRESENT**

BWX | 6-05-2025



No New Lesions!

Case Study: Ava N.
Pediatric Case



11-21-2019

FidaLab Caries | 11-2019

Score: 4.36, Moderate Risk

Results								
Type	Bacteria	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰
I	Streptococcus mutans	0.1	*					
	Streptococcus sobrinus	N/A						
II	Lactobacillus casei	16.2						
	Lactobacillus fermentum	9.4						
	Lactobacillus rhamnosus	3.5						
	Bifidobacterium dentium	285.0						
	Slackia exigua	N/A						
III	Scardovia wiggsiae	N/A						
	Propionibacterium acidifaciens	N/A						
	Actinomyces viscosus	61.1						
	Veillonella spp	16341.2						

N/A - Below detection limit * - <10⁵ copies / sample

Cariogenic Pathogen Composite Score: 4.36

Composite Score Key

<1 = Low Risk

> 1 but < 10 = Moderate Risk


> 10 = High Risk



April 2020 – November 2022: Cavity Management Planning Consultation

- Reviewed importance of fluoride and frequent radiographs for monitoring lesions
- Recommended homecare: **CariFree Treatment Rinse, CariFree CTx4 gel toothpaste, some nutritional counseling**
- Balance* by Dr. Kim Kutsch dispensed to parents for reading.
- All treatment & homecare recommendations reviewed at a **Cavity Management Planning Consultative session** with her parent.
- Seen every 6 months per patient request; fluoride varnish at each visit, annual bitewings for monitoring.**

LOIS : HVA 4-14-20



Dentistry for Total Body Wellness

Personalized Cavity Management Planner

Personalized Summary

- Eliminate or reduce the bacteria that causes cavities**
 - Mechanical Disruption:
 - Waterpik or Oral Irrigator
 - Power Brush
 - Dental Floss
 - Chemical Disruption with a home rinse:
 - CTx4 Disinfectant Treatment Rinse
 - Antibiotic Rinse
 - Bleach Based Home Remedy Rinse
 - Fluoride Rinse- ACT, etc.
 - Xylitol Rinse
 - Xylitol Products:
 - Xylitol Gum or Xylitol Lollipops
 - Xylitol Mints or Xylimelts
 - Powdered Xylitol in a Delivery System
 - Medicated Pastes or Gels:
 - Clinpro 5000 Paste
 - CTx4 Gel 5000
 - Daily Defense Paste
 - CTx4 1100 Gel
 - CTx3 Gel
 - Tray Carrier to Increase Contact time
 - Other Paste/Gel
- Eliminate sugar and use a substitute bacteria can't use**
 - Read the books *SOS, Balance, and Fat Chance*
 - Switch to Xylitol as a sweetener *Erythritol*
 - Avoid hidden sugar sources like- Tums, Juices, Cough Drops, Medications, White Bread, Peanut Butter, School Lunches, etc.
- Use products that reduce mouth acidity and restore natural pH**
 - Consider a mouthwash with a pH above 7.0 - *Sheet given*
 - Buy alkaline water such as Essentia or *Ice Mountain*
 - Switch to Xylitol as a sweetener
 - Neutralize mouth after meals with Xylitol or Alkaline products
- Add foods to your diet that don't allow bacteria to grow**
 - Eliminate foods containing refined sugar especially high fructose corn syrup
 - Eliminate juices with added sugar and soft drinks
 - Eat foods with alkaline properties- green leafy vegetables, foods high in arginine
- Disrupt the bacteria and harden teeth with the proper paste**
 - CariFree CTx4 gel
 - Clinpro 5000 Paste
 - CariFree CTx3 gel
 - Daily Defense Paste
 - Other Paste/Gel

Final thoughts

recommend 3x week - high fluoride paste CTx4 gel 5000

© 2017 Wellness Dentistry Network
Douglas Thompson, D.D.S 3684 W. Maple Rd. Bloomfield Hills, MI 48301 248-642-1000 www.ioralmed.com

FidaLab Caries | 08-2020

Score: .33, Low Risk

Results								
Type	Bacteria	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰
I	Streptococcus mutans	N/A						
	Streptococcus sobrinus	N/A						
II	Lactobacillus casei	N/A						
	Lactobacillus fermentum	0.0 *						
	Lactobacillus rhamnosus	3.1						
	Bifidobacterium dentium	N/A						
	Slackia exigua	N/A						
III	Scardovia wiggisiae	N/A						
	Propionibacterium acidifaciens	N/A						
	Actinomyces viscosus	31.3						
	Veillonella spp	14826.7						

N/A - Below detection limit

* - <10⁵ copies / sample






Cariogenic Pathogen Composite Score: 0.33

Composite Score Key

<1 = Low Risk

> 1 but < 10 = Moderate Risk

> 10 = High Risk






11-05-2020

10-27-2021

5-17-2022

Caries Case Study | Ava N. | Salivary Diagnostic Test Review



Ava's FIDA Lab results after three years of maintenance

Composite Score Key

- <1 = Low Risk
- > 1 but < 10 = Moderate Risk
- > 10 = High Risk





November 2022

- Ava's last appointment in our office was in **November 2022 due to her family moving 1 hour away from our office.**
- She planned to continue dental care at a clinic closer to her new home.

January 2025

- In early 2025, we contacted Ava's parents to follow up on her care.
- **By January 2025 (age 14), we completed a new Fidalab Comprehensive Caries test, bitewing radiographs, and clinical photos.**





Caries Case Study | Ava N. | Salivary Diagnostic 5-Year Test Review

			Results	
Type	Bacteria	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶
I	Streptococcus mutans	0.1*		
	Streptococcus sobrinus	N/A		
	Lactobacillus casei	16.2		
II	Lactobacillus fermentum	9.4		
	Lactobacillus rhamnosus	3.5		
	Bifidobacterium dentium	285.0		
	Slackia exigua	N/A		
III	Scardovia wiggisiae	N/A		
	Propionibacterium acidifaciens	N/A		
	Actinomyces viscosus	61.1		
	Veillonella spp	16341.2		

11-2019
Score: 4.36, Moderate

N/A - Below detection limit * - <10⁵ copies / sample Cariogenic Pathogen Composite Score: 4.36

			Results	
Type	Bacteria	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶
I	Streptococcus mutans	N/A		
	Streptococcus sobrinus	N/A		
	Lactobacillus casei	N/A		
II	Lactobacillus fermentum	0.0*		
	Lactobacillus rhamnosus	3.1		
	Bifidobacterium dentium	N/A		
	Slackia exigua	N/A		
III	Scardovia wiggisiae	N/A		
	Propionibacterium acidifaciens	N/A		
	Actinomyces viscosus	31.3		
	Veillonella spp	14826.7		

08-2020
Score: .33, Low

N/A - Below detection limit * - <10⁵ copies / sample Cariogenic Pathogen Composite Score: 0.33

			Results	
Type	Organism	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶
A	Streptococcus mutans	N/A		
	Streptococcus sobrinus	N/A		
B	Lactobacillus casei	N/A		
	Lactobacillus fermentum	00.010	*	
	Lactobacillus rhamnosus	N/A		
	Bifidobacterium dentium	00.474	*	
C	Slackia exigua	N/A		
	Actinomyces viscosus	310.344		
	Veillonella spp	14263.810		
	Scardovia wiggisiae	N/A		
	Propionibacterium acidifaciens	N/A		

11-2022
Score .11, Low

N/A Below Detection Limit * - < 10⁵ Copies/sample Cariogenic Pathogen Composite Score:00.11

FidaLab Comprehensive Caries | 01-2025
Score: .02 Low Risk | No Yeast Detected

A. Bacteria Pathogens Results								
Type	Organism	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰
A	Streptococcus mutans	N/A						
	Streptococcus sobrinus	N/A						
B	Lactobacillus casei	N/A						
	Lactobacillus fermentum	N/A						
	Lactobacillus rhamnosus	N/A						
	Bifidobacterium dentium	N/A						
C	Slackia exigua	N/A						
	Actinomyces viscosus	00.960	*					
	Veillonella spp	2498.372						
	Scardovia wiggisiae	N/A						
	Propionibacterium acidifaciens	N/A						

N/A Below Detection Limit * - < 10⁵ Copies/sample Bacterial pathogens composite score:00.02

B. Fungal Pathogens Results	
Candida albicans	Not Detected
Candida dubliniensis	Not Detected
Candida parapsilosis	Not Detected
Candida glabrata	Not Detected
Candida krusei	Not Detected
Candida tropicalis	Not Detected
Candida kefyr	Not Detected
Candida guilliermondii	Not Detected
Candida lusitanae	Not Detected
Candida rugosa	Not Detected

			Results	
Type	Bacteria	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶
I	Streptococcus mutans	0.1*		
	Streptococcus sobrinus	N/A		
	Lactobacillus casei	16.2		
II	Lactobacillus fermentum	9.4		
	Lactobacillus rhamnosus	3.5		
	Bifidobacterium dentium	285.0		
	Slackia exigua	N/A		
III	Scardovia wiggisiae	N/A		
	Propionibacterium acidifaciens	N/A		
	Actinomyces viscosus	61.1		
	Veillonella spp	16341.2		

11-2019
Score: 4.36, Moderate

N/A - Below detection limit * - <10⁵ copies / sample

Cariogenic Pathogen Composite Score: 4.36

			Results	
Type	Bacteria	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶
I	Streptococcus mutans	N/A		
	Streptococcus sobrinus	N/A		
	Lactobacillus casei	N/A		
II	Lactobacillus fermentum	0.0*		
	Lactobacillus rhamnosus	3.1		
	Bifidobacterium dentium	N/A		
	Slackia exigua	N/A		
III	Scardovia wiggisiae	N/A		
	Propionibacterium acidifaciens	N/A		
	Actinomyces viscosus	31.3		
	Veillonella spp	14826.7		

08-2020
Score: .33, Low

N/A - Below detection limit * - <10⁵ copies / sample

Cariogenic Pathogen Composite Score: 0.33

			Results	
Type	Organism	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶
A	Streptococcus mutans	N/A		
	Streptococcus sobrinus	N/A		
B	Lactobacillus casei	N/A		
	Lactobacillus fermentum	00.010	*	
	Lactobacillus rhamnosus	N/A		
	Bifidobacterium dentium	00.474	*	
C	Slackia exigua	N/A		
	Actinomyces viscosus	310.344		
	Veillonella spp	14263.810		
	Scardovia wiggisiae	N/A		
	Propionibacterium acidifaciens	N/A		

11-2022
Score .11, Low

N/A Below Detection Limit * - < 10⁵ Copies/sample

Cariogenic Pathogen Composite Score:00.11

FidaLab Comprehensive Caries | 01-2025
Score: .02 Low Risk | No Yeast Detected

A. Bacteria Pathogens Results								
Type	Organism	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰
A	Streptococcus mutans	N/A						
	Streptococcus sobrinus	N/A						
B	Lactobacillus casei	N/A						
	Lactobacillus fermentum	N/A						
	Lactobacillus rhamnosus	N/A						
	Bifidobacterium dentium	N/A						
C	Slackia exigua	N/A						
	Actinomyces viscosus	00.960	*					
	Veillonella spp	2498.372						
	Scardovia wiggisiae	N/A						
	Propionibacterium acidifaciens	N/A						

N/A Below Detection Limit * - < 10⁵ Copies/sample

Bacterial pathogens composite score:00.02

B. Fungal Pathogens Results	
<i>Candida albicans</i>	Not Detected
<i>Candida dubliniensis</i>	Not Detected
<i>Candida parapsilosis</i>	Not Detected
<i>Candida glabrata</i>	Not Detected
<i>Candida krusei</i>	Not Detected
<i>Candida tropicalis</i>	Not Detected
<i>Candida kefyr</i>	Not Detected
<i>Candida guilliermondii</i>	Not Detected
<i>Candida lusitanae</i>	Not Detected
<i>Candida rugosa</i>	Not Detected

No New Lesions!


Case Study: Daniel M.
Traditional Restorative Method

Medical / Dental History Review

- Initially presented in our office in 2018
- Estimated health as 'Excellent'
- **Had not been to a dentist in 20 years**
- **Reported pain in lower left; patient suspected a cavity.**



FMX | 03-2018



No obvious pathology noted. Diagnosis is debris around the 3rd molar on left lower side causing pain. Pain was resolved with irrigation.

Case Narrative

- Patient then did not return to our office until July 2025.



July 2025

- July 2025: Patient presented for a cleaning.
 - **Carious lesions present on almost all permanent teeth with heavy staining. Erosion present on occlusal surfaces.**
 - Patient reports no changes to his health, medications, home habits, alcohol/drug use, or diet from previously stated in 2018.
 - Recommended a biofilm assessment.



CariScreen & Caries Risk Assessment, 7-15-2025

Score: 7777, High Risk

CRA Form Adults and Children 6+

Four-Step Guide to Cavity Prevention

Name: Daniel Date: 7/15/25

1 Risk & Disease Discovery Complete the section below by circling the most appropriate answer for each question.

RISK FACTORS		
SALIVA		
Do you take medications daily? If so, how many?	NO	YES (5)
Do you feel as though you have a dry mouth at any time of the day or night?	NO	YES
DIET		
Do you drink liquids other than water more than 2 times daily between meals?	NO	YES
Do you snack daily between meals?	NO	YES
BIOFILM		
Do you notice plaque build-up on your teeth between brushings?	NO	YES
CariScreen reading results:	LOW (0-1500)	HIGH (1501-9999) 7777
Disease Indicators		
New/Progressing visible cavitations?	NO	YES
New/Progressing approximal radiographic radiolucencies?	NO	YES
New/Active white spot lesions?	NO	YES
Is decay history a concern?	NO	YES

Healthy	* Risk Factors	* Disease Indicators
1 - Low Risk	2 - Moderate Risk	3 - High Risk
CDT Code D0601	CDT Code D0602	CDT Code D0603

Once risk level diagnosis has been made, refer to next page for treatment options and reassessment protocol.

CARIFREE® | CRM®

CRA Form Adults and Children 6+

Four-Step Guide to Cavity Prevention

Name: Daniel Date: 7/15/25

3 Treatment Options

1 - Low Risk	2 - Moderate Risk	3 - High Risk

Dry Mouth
If dry mouth diagnosis has been made or if answered yes to dry mouth question, it is recommended to add moisturizing products such as Spray, Breath Mints, or Fluoride Free Gel for use throughout the day.

4 Reassessment Protocol
High Risk: Reassess every three months until moderate risk.
Moderate Risk: Reassess at every re-care appointment until low risk.
Low Risk: Reassess annually.

Four-Step Guide to Cavity Prevention

CARIFREE® | CRM®

CRA FORM Adults and Children Age 6+

First name: Daniel Last name: Date: 7/15/25

Due to new research on cavities and what causes them, we know everyone is at risk of developing decay at some point during their lifetime. The goal of this assessment form and the bacterial screening test is to determine your likelihood of experiencing new decay in the next 12 months. Please fill out the "Patient Use" section of this form to the best of your ability. These items will be discussed with your dental professional during your appointment today. Questions about this form? See the back for Q&A.

Would you like a free bacterial screening test to help determine your risk for cavities? (The test is a quick, painless swab of your teeth.)

If diagnosed at risk for cavities today, would you be interested in discussing treatment options?

If needed, are you willing to modify your dietary habits?

	yes	maybe	no
Would you like a free bacterial screening test to help determine your risk for cavities?	yes		no
If diagnosed at risk for cavities today, would you be interested in discussing treatment options?	yes	maybe	no
If needed, are you willing to modify your dietary habits?	yes	maybe	no

RISK FACTORS

	no	yes
Do you notice plaque build-up on your teeth between brushings?	no	yes
Do you take medications daily? If yes, how many? (#)	no	yes
Do you feel like you have a dry mouth at any time of the day or night?	no	yes
Do you drink liquids other than water more than 2 times daily between meals?	no	yes
Do you snack daily between meals?	no	yes
Do you have oral appliances present?	no	yes
Do any of these other health concerns apply to you? (check all that apply)		yes
<input type="checkbox"/> Frequent tobacco use	no	yes
<input type="checkbox"/> Acid reflux		
<input type="checkbox"/> Diabetes		
<input type="checkbox"/> Head/neck radiation therapy		
<input type="checkbox"/> Other drug use		
<input type="checkbox"/> Bulimia		
<input type="checkbox"/> Sjogren's Syndrome		

DISEASE INDICATORS

	no	yes
New/Progressing Visible Cavitations	no	yes
New/Progressing Approximal Radiographic Radiolucencies	no	yes
New/Active White Spot Lesions	no	yes
Decay History is a Concern	no	yes

BIOFILM CHALLENGE
CariScreen Bacterial Assessment (0-1500 low, 1501-9999 high)

	low	high
CariScreen Bacterial Assessment (0-1500 low, 1501-9999 high)	low	high 7777

PROFESSIONAL ASSESSMENT SUMMARY

	no	yes
Risk Factors are a Concern	no	yes
Disease Indicators are a Concern	no	yes
Biofilm Challenge is a Concern	no	yes

RISK IDENTIFICATION Transfer information above to boxes below to determine risk.

	LOW RISK	MODERATE RISK	HIGH RISK	HIGH RISK	HIGH/EXTREME RISK
<input type="checkbox"/> Risk Factors					
<input type="checkbox"/> Disease Indicators					
<input type="checkbox"/> Biofilm Challenge					
<input type="checkbox"/> Risk Factors					
<input type="checkbox"/> Disease Indicators					
<input type="checkbox"/> Biofilm Challenge					
<input type="checkbox"/> Risk Factors					
<input type="checkbox"/> Disease Indicators					
<input type="checkbox"/> Biofilm Challenge					
<input type="checkbox"/> Risk Factors					
<input type="checkbox"/> Disease Indicators					
<input type="checkbox"/> Biofilm Challenge					

1 2 3 4 5

RECOMMENDED | PROVISIONAL | DECLINE

CARIFREE® | CRM®

FidaLab Comprehensive Caries, Pre-Restorative | 07-2025

Score: 311.10, High Risk | **No Yeast Detected**

A. Bacteria Pathogens Results								
Type	Organism	Copies/oral rinse sample (x10 ^{^5})	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰
A	Streptococcus mutans	01.530						
	Streptococcus sobrinus	N/A						
B	Lactobacillus casei	3416.959						
	Lactobacillus fermentum	07.460						
	Lactobacillus rhamnosus	823.236						
	Bifidobacterium dentium	134.523						
	Slackia exigua	N/A						
C	Actinomyces viscosus	13.143						
	Veillonella spp	11829.211						
	Scardovia wiggisiae	13.592						
	Propionibacterium acidifaciens	19.231						

N/A Below Detection Limit

* < 10⁵ Copies/sample

Bacterial pathogens composite score:311.10

B. Fungal Pathogens Results	
<i>Candida albicans</i>	Not Detected
<i>Candida dubliniensis</i>	Not Detected
<i>Candida parapsilosis</i>	Not Detected
<i>Candida glabrata</i>	Not Detected
<i>Candida krusei</i>	Not Detected
<i>Candida tropicalis</i>	Not Detected
<i>Candida kefyr</i>	Not Detected
<i>Candida guilliermondii</i>	Not Detected
<i>Candida lusitanae</i>	Not Detected
<i>Candida rugosa</i>	Not Detected

Saliva-Check Buffer Test | 07-2025

Saliva-Check Buffer Tests:

1. Resting saliva → **Normal (>60 seconds)**
 - Visual inspection of level of hydration; observe lower lip labial mucosa gland secretion.
2. Saliva consistency → **Normal (Watery clear saliva)**
 - Visually assessing resting salivary consistency.
3. pH measurement → **Moderately acidic (6.0 – 6.6 pH)**
 - pH testing strip into salivary sample for 10 seconds
4. Saliva quantity → **Normal (>5.0 mL)**
 - Stimulate saliva for 30 seconds by chewing piece of wax. Expectorate saliva into cup and continue chewing for 5 minutes with expectoration.
5. Buffering capacity: **Very low (0 – 5 score)**
 - Pipette saliva from collection cup onto 3 test pads, wait 2 minutes, and read test pad.



Caries Case Study | Daniel M. | Initial Salivary Diagnostic Test, Results Review

A. Bacteria Pathogens Results								
Type	Organism	Copies/oral rinse sample (x10 ⁵)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰
A	Streptococcus mutans	01.530	█					
	Streptococcus sobrinus	N/A		█				
B	Lactobacillus casei	3416.959	█					
	Lactobacillus fermentum	07.460	█					
	Lactobacillus rhamnosus	823.236	█					
	Bifidobacterium dentium	134.523	█					
	Slackia exigua	N/A						
C	Actinomyces viscosus	13.143	█					
	Veillonella spp	11829.211	█					
	Scardovia wiggsiae	13.592	█					
	Propionibacterium acidifaciens	19.231	█					

N/A Below Detection Limit

* < 10⁵ Copies/sample

Bacterial pathogens composite score:311.10

TEST RESULTS SUMMARY

1. Dysbiotic microbiome with high levels of cariogenic bacteria present
2. Dan's saliva is naturally more acidic; 6.0 – 6.6 pH
3. Dan's saliva naturally has very low buffering capabilities

CRA Form Adults and Children 6- Four-Step Guide to Cavity Prevention

Name: Daniel Date: 7/15/25

1 Risk & Disease Discovery Complete the section below by circling the most appropriate answer for each question.

RISK FACTORS	
SALIVA	
Do you take medications daily? If so, how many?	NO YES (5)
Do you feel as though you have a dry mouth at any time of the day or night?	NO YES
DIET	
Do you drink liquids other than water more than 2 times daily between meals?	NO YES
Do you snack daily between meals?	NO YES
BIOFILM	
Do you notice plaque build-up on your teeth between brushings?	NO YES
CarScreen reading results:	LOW (0-1000) HIGH (101-9999)
DISEASE INDICATORS	
New/Progressing visible cavitations?	NO YES
New/Progressing approximal radiographic radioluscencies?	NO YES
New/Active white spot lesions?	NO YES
Is decay history a concern?	NO YES

2 Risk Identification Determine risk based on answers above; one blue response indicates moderate risk, one red response indicates high risk.

Healthy	+ Risk Factors	+ Disease Indicators
1 - Low Risk	2 - Moderate Risk	3 - High Risk
CDT Code D0601	CDT Code D0602	CDT Code D0603

Once risk level diagnosis has been made, refer to next page for treatment options and reassessment protocol.

CARIFREE® | CRM® v16



Case Narrative: August 2025

- All current restorative needs were addressed in July
- Only homecare recommendations followed during restorative treatment was utilizing **CariFree ProGel 5000**
- Did not start any biofilm management yet. Wanted to do restorations and reevaluate the microbiome.
- **Fidalab Comprehensive Caries test re-administered**



Caries Case Study | Daniel M. | Salivary Diagnostic Test, Pre vs Post Restorative

FidaLab Comprehensive Caries, Pre-Restorative | 07-2025

Score: 311.10, High Risk | **No Yeast Detected**

A. Bacteria Pathogens Results									
Type	Organism	Copies/oral rinse sample (x10 ⁴ %)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰	10 ¹¹
A	Streptococcus mutans	01.530							
	Streptococcus sobrinus	N/A							
B	Lactobacillus casei	3416.959							
	Lactobacillus fermentum	07.460							
	Lactobacillus rhamnosus	823.236							
	Bifidobacterium dentium	134.523							
	Slackia exigua	N/A							
C	Actinomyces viscosus	13.143							
	Veillonella spp	11829.211							
	Scardovia wiggisiae	13.592							
	Propionibacterium acidifaciens	19.231							

N/A Below Detection Limit * < 10⁵ Copies/sample

Bacterial pathogens composite score:311.10

B. Fungal Pathogens Results	
Candida albicans	Not Detected
Candida dubliniensis	Not Detected
Candida parapsilosis	Not Detected
Candida glabrata	Not Detected
Candida krusei	Not Detected
Candida tropicalis	Not Detected
Candida kefir	Not Detected
Candida guilliermondii	Not Detected
Candida lusitanae	Not Detected
Candida rugosa	Not Detected

FidaLab Comprehensive Caries, Post-Restorative | 08-2025

Score: 373.95, High Risk | **No Yeast Detected**

A. Bacteria Pathogens Results									
Type	Organism	Copies/oral rinse sample (x10 ⁴ %)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰	10 ¹¹
A	Streptococcus mutans	01.100							
	Streptococcus sobrinus	N/A							
B	Lactobacillus casei	1577.486							
	Lactobacillus fermentum	356.552							
	Lactobacillus rhamnosus	2891.586							
	Bifidobacterium dentium	3038.298							
	Slackia exigua	N/A							
C	Actinomyces viscosus	986.899							
	Veillonella spp	12283.857							
	Scardovia wiggisiae	07.421							
	Propionibacterium acidifaciens	19.026							

Bacterial pathogens composite score:373.95

B. Fungal Pathogens Results	
Candida albicans	Not Detected
Candida dubliniensis	Not Detected
Candida parapsilosis	Not Detected
Candida glabrata	Not Detected
Candida krusei	Not Detected
Candida tropicalis	Not Detected
Candida kefir	Not Detected
Candida guilliermondii	Not Detected
Candida lusitanae	Not Detected
Candida rugosa	Not Detected

HIGHER SCORE?

Caries Case Study | Daniel M. | Salivary Diagnostic Test, Pre vs Post Restorative

FidaLab Comprehensive Caries, Pre-Restorative | 07-2025

Score: 311.10, High Risk | **No Yeast Detected**

A. Bacteria Pathogens Results

Type	Organism	Copies/oral rinse sample (x10 ⁴ %)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰
A	Streptococcus mutans	01.530	█					
	Streptococcus sobrinus	N/A						
B	Lactobacillus casei	3416.959	█	█	█	█		
	Lactobacillus fermentum	07.460	█					
	Lactobacillus rhamnosus	823.236	█	█	█			
	Bifidobacterium dentium	134.523	█	█	█			
	Slackia exigua	N/A						
C	Actinomyces viscosus	13.143	█					
	Veillonella spp	11829.211	█	█	█	█		
	Scardovia wiggsiae	13.592	█					
	Propionibacterium acidifaciens	19.231	█					

N/A Below Detection Limit * < 10⁵ Copies/sample

Bacterial pathogens composite score:311.10

B. Fungal Pathogens Results

<i>Candida albicans</i>	Not Detected
<i>Candida dubliniensis</i>	Not Detected
<i>Candida parapsilosis</i>	Not Detected
<i>Candida glabrata</i>	Not Detected
<i>Candida krusei</i>	Not Detected
<i>Candida tropicalis</i>	Not Detected
<i>Candida kefyri</i>	Not Detected
<i>Candida guilliermondii</i>	Not Detected
<i>Candida lusitanae</i>	Not Detected
<i>Candida rugosa</i>	Not Detected

FidaLab Comprehensive Caries, Post-Restorative | 08-2025

Score: 373.95, High Risk | **No Yeast Detected**

A. Bacteria Pathogens Results

Type	Organism	Copies/oral rinse sample (x10 ⁴ %)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰
A	Streptococcus mutans	01.100	█					
	Streptococcus sobrinus	N/A						
B	Lactobacillus casei	1577.486	█	█	█	█		
	Lactobacillus fermentum	356.552	█	█	█			
	Lactobacillus rhamnosus	2891.586	█	█	█	█		
	Bifidobacterium dentium	3038.298	█	█	█	█		
	Slackia exigua	N/A						
C	Actinomyces viscosus	986.899	█	█	█			
	Veillonella spp	12283.857	█	█	█	█		
	Scardovia wiggsiae	07.421	█					
	Propionibacterium acidifaciens	19.026	█					

Bacterial pathogens composite score:373.95

B. Fungal Pathogens Results

<i>Candida albicans</i>	Not Detected
<i>Candida dubliniensis</i>	Not Detected
<i>Candida parapsilosis</i>	Not Detected
<i>Candida glabrata</i>	Not Detected
<i>Candida krusei</i>	Not Detected
<i>Candida tropicalis</i>	Not Detected
<i>Candida kefyri</i>	Not Detected
<i>Candida guilliermondii</i>	Not Detected
<i>Candida lusitanae</i>	Not Detected
<i>Candida rugosa</i>	Not Detected

HIGHER SCORE?

WHAT DOES THIS MEAN?

The traditional dental strategy of performing restorative work and prescribing supplemental fluoride fails to effectively target one of the main etiologies of dental caries, which is the cariogenic biofilm.

Case Study: Peter R.
Wellness Dentistry Method

Medical / Dental History Review

- 66-year-old male.
- Reports overall health as "good."
- Light sleeper.
- No history of sleep apnea.
- Borderline high cholesterol.
- **Average arterial buildup; wants to reverse heart disease.**
- Takes dietary supplements for heart health
- **Congenitally missing teeth #7 and #10.**
- **Referred by physician and reports frequent cavities**
- Issues with root canal on tooth #15.
- Had completed CIMT tests, OralDNA MyPerioPath, previous sleep test, pulmonary function reports, and bloodwork sent to our office.

What is your estimate of your general health?

Excellent	No
Good	YES
Fair	No
Poor	No

14. Chronic ear infections, tuberculosis, measles, chicken pox	YES	Not sure, but I may have had measles and or chicken pox when I was a child. Did they have vaccines in 60's?
15. Breathing problems (e.g. asthma, stuffy nose, sinus congestion)	No	
16. Sleep problems (e.g. sleep apnea, snoring, insomnia, restless sleep, bedwetting)	YES	I am a very light sleeper, wake frequently majority of the time, but do not have sleep apnea (was tested)
17. Kidney disease	No	
18. Liver disease or jaundice	No	
19. Vertigo (e.g. "the room is spinning")	No	
20. Thyroid, parathyroid disease, or calcium deficiency	No	
21. Hormone deficiency or imbalance (e.g. polycystic ovarian syndrome)	No	
22. High cholesterol or taking statin drugs	YES	Cholesterol is on borderline for high
48. Presently being treated for any other illness	YES	Trying to reduce heart disease, I am average as far as the level of arterial buildup, but trying to reverse/stop it from progressing.
49. Aware of any change in your health in the last 24 hours (e.g. fever, chills, new cough, or diarrhea)	No	
50. Taking medication for weight management	No	
51. Taking dietary supplements, vitamins, and/or probiotics	YES	Methyl B Complex, K-Force (Vit k and D), CoQ10, EPA, Magnesium Glycinate.
52. Often exhausted or fatigued	No	

FMX | 03-2025



**Treatment plan: #2 crown, #3-MO, #15 crown due to recurrent decay on root surface, #18 RCT and crown.
Awaiting Fidalab Comprehensive Caries results for further homecare recommendations**

FidaLab Comprehensive Caries | 03-2025

Score: 301.12, High Risk | *C. albicans* Detected

A. Bacteria Pathogens Results								
Type	Organism	Copies/oral rinse sample (x10 ⁶)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰
A	<i>Streptococcus mutans</i>	52.413						
	<i>Streptococcus sobrinus</i>	N/A						
B	<i>Lactobacillus casei</i>	737.884						
	<i>Lactobacillus fermentum</i>	1435.154						
	<i>Lactobacillus rhamnosus</i>	1335.157						
	<i>Bifidobacterium dentium</i>	1038.053						
	<i>Slackia exigua</i>	N/A						
C	<i>Actinomyces viscosus</i>	4929.965						
	<i>Veillonella spp</i>	37408.757						
	<i>Scardovia wiggisiae</i>	74.059						
	<i>Propionibacterium acidifaciens</i>	N/A						

N/A Below Detection Limit * -< 10⁵ Copies/sample

Bacterial pathogens composite score:301.12

<i>Candida albicans</i>	Detected
<i>Candida dubliniensis</i>	Not Detected
<i>Candida parapsilosis</i>	Not Detected
<i>Candida glabrata</i>	Not Detected
<i>Candida krusei</i>	Not Detected
<i>Candida tropicalis</i>	Not Detected
<i>Candida kefyr</i>	Not Detected
<i>Candida guilliermondii</i>	Not Detected
<i>Candida lusitanae</i>	Not Detected
<i>Candida rugosa</i>	Not Detected

Case Narrative, Cont'd

- 03/2025: Completed #2 crown and #3 MO restoration.
- 04/2025: Patient seen for caries consultation at crown seat appointment.
 - **High caries risk due to composite score of 301.12, presence of *S. mutans* and *C. albicans***
- **Prescribed a 1000 mg/tsp amoxicillin/nystatin rinse to be used for 2 weeks to manage dysbiotic oral environment. Follow with CariFree Treatment Rinse twice daily for one minute.**
- **Follow-up caries test scheduled for 6 weeks post-rinse.**
- Ordered **CHL lab panel** for blood analysis.
- Patient purchased **CariFree Treatment Rinse and a Waterpik**. Dispensed & personalized **Order of Operations** take-home resource



Recommended Home Care Order of Operations

Mastering the correct order of your homecare routine lays the foundation for great oral care. It ensures that each step, whether it's flossing, brushing, or using mouthwash, is strategically placed to achieve maximum efficacy. By adhering to this sequence, you not only optimize the effectiveness of your oral hygiene products but also enhances the overall health of your teeth and gums, promoting a confident and radiant smile.

Optimize Your Personalized Homecare Routine

For optimal clinical outcomes, follow this precise order in your homecare routine:

1. **In-Between Tooth Care:** *Waterpik*
 - Start with either using a Waterpik, floss, or any "in between" the teeth cleaning method that is effective to remove debris between the teeth. *(Using 2-5 mL of bleach inside reservoir)*
2. **Mouthrinse:** *Amox/Nystatin (2wks), then CariFree rinse.*
 - If recommended, rinse your mouth with desired rinse for **one minute**. If you have a healthy balanced mouth, you should not need a mouthrinse on a regular basis. If your goal is to freshen your breath, consider xylitol-based gum or mints.
3. **Brush:** *W/ CTx4 gel purchased @ last visit*
 - Brush all surfaces of your teeth with your prescribed toothpaste for a **full two minutes**.
4. **Spit, Don't Rinse:**
 - After brushing, spit out excess toothpaste without rinsing. This maximizes the absorption of the active ingredients in the toothpaste.

Consistency in your homecare routine is a key determinant of long-term oral health outcomes. Small, daily efforts—when performed correctly—can create impactful changes in the oral environment, influencing everything from microbial balance to gum health and cavity control. Your adherence to these steps plays a crucial role in maintaining a stable, disease-resistant oral cavity!

© 2025 Wellness Dentistry Network

06-2025

A. Bacteria Pathogens Results

Type	Organism	Copies/oral rinse sample (x10 ⁶)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰
A	Streptococcus mutans	01.833						
	Streptococcus sobrinus	N/A						
B	Lactobacillus casei	N/A						
	Lactobacillus fermentum	N/A						
	Lactobacillus rhamnosus	N/A						
	Bifidobacterium dentium	01.076						
	Slackia exigua	00.046	*					
C	Actinomyces viscosus	224.474						
	Veillonella spp	31210.178						
	Scardovia wiggsiae	27.009						
	Propionibacterium acidifaciens	N/A						

N/A Below Detection Limit

* < 10⁵ Copies/sample

Bacterial pathogens composite score:01.57

B. Fungal Pathogens Results

<i>Candida albicans</i>	Not Detected
<i>Candida dubliniensis</i>	Not Detected
<i>Candida parapsilosis</i>	Not Detected
<i>Candida glabrata</i>	Not Detected
<i>Candida krusei</i>	Not Detected
<i>Candida tropicalis</i>	Not Detected
<i>Candida kefyr</i>	Not Detected
<i>Candida guilliermondii</i>	Not Detected
<i>Candida lusitanae</i>	Not Detected
<i>Candida rugosa</i>	Not Detected

FidaLab Comprehensive Caries | 03-2025
Score: 301.12, High Risk |
C. albicans Detected

FidaLab Comprehensive Caries | 06-2025
Score: 1.57, Low Risk
No Yeast Detected

10-2025

A. Bacteria Pathogens Results

Type	Organism	Copies/oral rinse sample (x10 ⁶)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰
A	Streptococcus mutans	00.596	*					
	Streptococcus sobrinus	N/A						
B	Lactobacillus casei	N/A						
	Lactobacillus fermentum	N/A						
	Lactobacillus rhamnosus	00.130	*					
	Bifidobacterium dentium	00.792	*					
	Slackia exigua	N/A						
C	Actinomyces viscosus	104.748						
	Veillonella spp	17003.961						
	Scardovia wiggisiae	28.965						
	Propionibacterium acidifaciens	N/A						

N/A Below Detection Limit

* < 10⁵ Copies/sample

Bacterial pathogens composite score:00.57

B. Fungal Pathogens Results

Candida albicans	Not Detected
Candida dubliniensis	Not Detected
Candida parapsilosis	Not Detected
Candida glabrata	Not Detected
Candida krusei	Not Detected
Candida tropicalis	Not Detected
Candida kefyr	Not Detected
Candida guilliermondii	Not Detected
Candida lusitanae	Not Detected
Candida rugosa	Not Detected

No New Lesions!

**FidaLab Comprehensive Caries
03-2025**

Score: 301.12, High Risk |
C. albicans Detected

**FidaLab Comprehensive Caries
06-2025**

Score: 1.57, Moderate Risk
No Yeast Detected

**FidaLab Comprehensive Caries
10-2025**

Score: 0.57, Low Risk
No Yeast Detected

10-2025

A. Bacteria Pathogens Results

Type	Organism	Copies/oral rinse sample (x10 ⁶)	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰
A	Streptococcus mutans	00.596	*					
	Streptococcus sobrinus	N/A						
B	Lactobacillus casei	N/A						
	Lactobacillus fermentum	N/A						
	Lactobacillus rhamnosus	00.130	*					
	Bifidobacterium dentium	00.792	*					
	Slackia exigua	N/A						
C	Actinomyces viscosus	104.748						
	Veillonella spp	17003.961						
	Scardovia wiggisiae	28.965						
	Propionibacterium acidifaciens	N/A						

N/A Below Detection Limit

* -< 10⁵ Copies/sample

Bacterial pathogens composite score:00.57

B. Fungal Pathogens Results

<i>Candida albicans</i>	Not Detected
<i>Candida dubliniensis</i>	Not Detected
<i>Candida parapsilosis</i>	Not Detected
<i>Candida glabrata</i>	Not Detected
<i>Candida krusei</i>	Not Detected
<i>Candida tropicalis</i>	Not Detected
<i>Candida kefyr</i>	Not Detected
<i>Candida guilliermondii</i>	Not Detected
<i>Candida lusitanae</i>	Not Detected
<i>Candida rugosa</i>	Not Detected

WHAT DOES THIS MEAN?

By addressing Peter's restorative needs, adjusting his homecare, and shifting one of the main underlying etiological agents of caries development (his dysbiotic microbiome), we were successfully able to alter his biofilm to healthier metrics, supporting long-term oral and systemic health.

Practical implementation of caries disease management in the era of CAMBRA: a guide for the general dentist

Douglas A. Young, DDS, EdD, MBA, MS ■ Ryan L. Quock, DDS ■ Susie Goolsby, DDS, MSHA
 Martha McComas, RDH, MS ■ Steven Powell, DDS ■ Gregory G. Zeller, DDS, MS
 Allen Wong, DDS, EdD ■ Brian B. Nový, DDS

Although the philosophy of caries management by risk assessment (CAMBRA) is widely accepted as an essential part of comprehensive disease management, its implementation in daily clinical practice has historically been hindered by a lack of concise guidance that bridges the gap between scholarly research and the logistics of patient care. With the aim of guiding comprehensive disease management for dental caries by integrating available evidence-based systems such as the American Dental Association Caries Classification System (ADA CCS) and caries risk assessment (CRA), this review presents clinical strategies in a functional format. Clinical caries management strategies that integrate the ADA CCS with CRA were compiled and evaluated by clinicians representing a variety of practice settings. Examples of clinical protocols were selected to illustrate the use of CRA at the patient level and ADA CCS lesion classification at the tooth level to comprehensively manage caries lesions in a high-carries-risk population. Comprehensive caries management combines the clinical use of CRA and caries lesion classification to guide clinical decisions and to support tooth preservation in high-risk populations. These strategies should be individualized and fulfill evidence-based clinical practice criteria with oral health as the primary outcome. Person-centered comprehensive caries management integrating CRA and classification of caries lesions has the potential to improve patient outcomes by addressing the underlying causes of caries disease. Documenting and measuring the processes at the tooth, patient, and practice levels are essential for generating outcomes that will drive future approaches aimed at improving patient health.

Received: February 16, 2025
 Revised: May 25, 2025
 Accepted: June 30, 2025

Keywords: caries classification system, caries management, caries risk assessment, caries risk management, clinical protocols, dental caries, oral health

Published with permission of the Academy of General Dentistry.
 © Copyright 2025 by the Academy of General Dentistry.
 All rights reserved. For printed and electronic reprints of this article for distribution, please contact scsreprints@sherdan.com.

The respect for comprehensive disease management as a whole, and caries management by risk assessment (CAMBRA) as a leading philosophy, has continued to increase among payers, educators, clinicians, organized oral healthcare, and the US Health Resources and Services Administration. Yet in an era of evidence-based clinical decision-making, the adoption of these concepts by more clinicians seems to be stalled by a lack of practical translational strategies that are urgently needed for the improvement of patient health. This is not to say that the underlying science is lacking or that research has failed to show any positive health outcomes from these approaches. Instead, the situation highlights the hurdles that make it difficult to evaluate and adopt the vast array of available information in the form of effective clinical protocols. Given that dental caries is one of the most common diseases treated in dentistry, it could be argued that every general dentist functions as a practical cariologist, requiring actionable information. This summary of clinical cariology based on the philosophy of CAMBRA is meant to provide guidance for general dentists to implement disease management concepts in daily practice and improve the oral health of the patients they treat.

Terminology

The term *caries* by itself is open to several interpretations, and *caries* is colloquially used to describe both the disease and the resulting symptomatic lesion; this may create confusing variations in diagnosis and treatment. Adoption of high-level consensus terminology, such as that described jointly by the International Association for Dental Research Cariology Research Group, the European Organization for Caries Research, and the International Caries Consensus Collaboration is recommended for more precise diagnosis, treatment, and communication.¹⁻³ For example, the term *dental caries* is used to describe the actual disease process, while *caries lesion* describes the clinical manifestation of the disease; the 2 terms are closely related yet distinct.

Caries management is defined as "...actions taken to interfere with mineral loss at all stages of the caries disease, including nonoperative and operative interventions/treatments."¹ The term *nonoperative* refers to treating the lesion without tooth structure removal by employing the following strategies: surface protection (including barriers or sealants), biofilm modification, remineralization via chemotherapeutics, saliva supplementation, and behavioral modification, including changes in diet and lifestyle.

	Low caries risk	Moderate caries risk	High caries risk	Extreme caries risk
Determining caries risk: The checked item farthest to the right determines overall caries risk	<input type="checkbox"/> No active caries lesion or radiolucency (cavitated or noncavitated) during past 24 mo <input type="checkbox"/> Healthy lifestyle habits, including diet with little to no consumption of simple sugars, infrequent snacking, and no drug or alcohol misuse <input type="checkbox"/> No visible plaque	<input type="checkbox"/> No active caries lesion or radiolucency (cavitated or noncavitated) during past 12 mo 1-2 of the following: <input type="checkbox"/> Unhealthy lifestyle habits, including occasional (≤ 2 times per day) between-meal snacks of simple sugars, or drug or alcohol misuse <input type="checkbox"/> Inadequate oral hygiene or visible plaque <input type="checkbox"/> Wearing dental or orthodontic appliances <input type="checkbox"/> Susceptible pits and fissures <input type="checkbox"/> Exposed root surfaces <input type="checkbox"/> Saliva-reducing factors (medications, radiation, or systemic diseases)	<input type="checkbox"/> Any active caries lesion or radiolucency (cavitated or noncavitated) during the past 12 mo <input type="checkbox"/> High bacterial load by measurement or observation (heavy plaque) 3 or more of the following: <input type="checkbox"/> Unhealthy lifestyle habits, including frequent (≥ 3 times per day) between-meal snacks of simple sugars or drug or alcohol misuse <input type="checkbox"/> Inadequate oral hygiene <input type="checkbox"/> Wearing dental or orthodontic appliances <input type="checkbox"/> Susceptible pits and fissures <input type="checkbox"/> Exposed root surfaces <input type="checkbox"/> Saliva-reducing factors (medications, radiation, or systemic diseases)	<input type="checkbox"/> High caries risk with hyposalivation by observation or measurement <input type="checkbox"/> High dependency on others for care
Circle overall caries risk	Low caries risk	Moderate caries risk	High caries risk	Extreme caries risk
Recall interval	12 mo	6 mo	3 mo	3 mo
Radiographs	24-36 mo	18-24 mo	6-18 mo	6 mo (until no new caries lesions)

*This new tabular caries risk assessment form, adapted from the caries management by risk assessment form published in 2013 and revised in 2019, automatically incorporates the instructions and guidelines of those protocols.^{1,2,8} The checked item farthest to the right determines the overall caries risk level as well as the recommended recall and radiographic intervals. This form, completed after the comprehensive clinical examination, is meant to be fast and simple; clinicians should be able to complete it within 1 to 2 minutes. The user starts in the last column on the right and works to the left. Studies have validated that the instructions and guidelines used in this table lead to selection of the correct caries risk levels no matter who fills out the form.²⁸⁻²⁷

CRA FORM

CRM Strategy

Risk level	Home care recommendations	Recall interval, mo	Radiographs, mo
Low	Effective toothbrushing and interdental cleaning, oral hygiene instruction, individualized diet modification (eg, reduced frequency and exposure to sugary snacks), and, OTC fluoride toothpaste	12	24-36
Moderate	<i>All of the above plus:</i> • Xylitol chewing gum or mints throughout the day • 0.05% sodium fluoride mouthrinse after meals, if possible • Alternative therapy: 5000 ppm fluoride toothpaste twice daily (instead of OTC fluoride toothpaste) and no mouthrinse	6	18-24
High	<i>All the above plus:</i> • 5000 ppm fluoride toothpaste (instead of OTC toothpaste), morning and night • Antibacterial ^P and pH neutralization strategies (eg, 0.2% hypochlorite rinse) before bedtime or topical office application of povidone iodine and fluoride varnish.	3	6-18
Extreme	<i>All the above plus:</i> • pH neutralization (neutralizing sprays, rinses, lozenges), xylitol gum, and/or xylitol lollipops throughout the day • Calcium and phosphate supplementation (eg, nano hydroxyapatite) throughout the day and right before bedtime	3	6 (until no new caries lesions)

Practical implementation of caries disease management in the era of CAMBRA: a guide for the general dentist. Young DA, Quock RL, Goolsby S, et al. *Gen Dent*. 2025;73(6):56-66. Agd.org/generaldentistry

Caries Risk Assessment Form

Determining caries risk: The checked item farthest to the right determines overall caries risk	Low caries risk	Moderate caries risk	High caries risk	Extreme caries risk
	<input type="checkbox"/> No active caries lesion or radiolucency (cavitated or noncavitated) during past 24 mo <input type="checkbox"/> Healthy lifestyle habits, including diet with little to no consumption of simple sugars, infrequent snacking, and no drug or alcohol misuse <input type="checkbox"/> No visible plaque	<input type="checkbox"/> No active caries lesion or radiolucency (cavitated or noncavitated) during past 12 mo <i>1-2 of the following:</i> <input type="checkbox"/> Unhealthy lifestyle habits, including occasional (≤ 2 times per day) between-meal snacks of simple sugars, or drug or alcohol misuse <input type="checkbox"/> Inadequate oral hygiene or visible plaque <input type="checkbox"/> Wearing dental or orthodontic appliances <input type="checkbox"/> Susceptible pits and fissures <input type="checkbox"/> Exposed root surfaces <input type="checkbox"/> Saliva-reducing factors (medications, radiation, or systemic diseases)	<input type="checkbox"/> Any active caries lesion or radiolucency (cavitated or noncavitated) during the past 12 mo <input type="checkbox"/> High bacterial load by measurement or observation (heavy plaque) <i>3 or more of the following:</i> <input type="checkbox"/> Unhealthy lifestyle habits, including frequent (≥ 3 times per day) between-meal snacks of simple sugars or drug or alcohol misuse <input type="checkbox"/> Inadequate oral hygiene <input type="checkbox"/> Wearing dental or orthodontic appliances <input type="checkbox"/> Susceptible pits and fissures <input type="checkbox"/> Exposed root surfaces <input type="checkbox"/> Saliva-reducing factors (medications, radiation, or systemic diseases)	<input type="checkbox"/> High caries risk with hyposalivation by observation or measurement <input type="checkbox"/> High dependency on others for care
Circle overall caries risk	Low caries risk	Moderate caries risk	High caries risk	Extreme caries risk
Recall interval	12 mo	6 mo	3 mo	3 mo
Radiographs	24-36 mo	18-24 mo	6-18 mo	6 mo (until no new caries lesions)

^aThis new tabular caries risk assessment form, adapted from the caries management by risk assessment form published in 2013 and revised in 2019, automatically incorporates the instructions and guidelines of those protocols.^{11,25} The checked item farthest to the right determines the overall caries risk level as well as the recommended recall and radiographic intervals. This form, completed after the comprehensive clinical examination, is meant to be fast and simple; clinicians should be able to complete it within 1 to 2 minutes. The user starts in the last column on the right and works to the left. Studies have validated that the instructions and guidelines used in this table lead to selection of the correct caries risk levels no matter who fills out the form.²⁵⁻²⁷

Agd.org/generaldentistry; Young DA, Quock RL, Goolsby S, et al. Practical implementation of caries disease management in the era of CAMBRA: a guide for the general dentist. *Gen Dent.* 2025;73(6):56-66.

Home Care - Caries Risk Management Strategy

Risk level	Home care recommendations	Recall interval, mo	Radiographs, mo
Low	Effective toothbrushing and interdental cleaning, oral hygiene instruction, individualized diet modification (eg, reduced frequency and exposure to sugary snacks), and, OTC fluoride toothpaste	12	24-36
Moderate	<i>All of the above plus:</i> <ul style="list-style-type: none"> • Xylitol chewing gum or mints throughout the day • 0.05% sodium fluoride mouthrinse after meals, if possible • Alternative therapy: 5000 ppm fluoride toothpaste twice daily (instead of OTC fluoride toothpaste) and no mouthrinse 	6	18-24
High	<i>All the above plus:</i> <ul style="list-style-type: none"> • 5000 ppm fluoride toothpaste (morning and night) • Antibacterial mouthrinse (chlorhexidine or povidone iodine) before bedtime or topical fluoride application or povidone iodine and fluoride varnish. 	3	6-18
Extreme	<i>All the above plus:</i> <ul style="list-style-type: none"> • pH neutralization (neutralizing toothpaste or xylitol lollipops throughout the day) • Calcium and phosphate supplements throughout the day and right before bedtime 	3	6 (until no new caries lesions)

BACTERIAL ASSESSMENT

SALIVA ASSESSMENT

Adapted from University of the Pacific AEGD clinics.

Agd.org/generaldentistry; Young DA, Quock RL, Goolsby S, et al. Practical implementation of caries disease management in the era of CAMBRA: a guide for the general dentist. *Gen Dent.* 2025;73(6):56-66.

ADA CCS

JADA Feb 2015

COVER STORY

The American Dental Association Caries Classification System for Clinical Practice

A report of the American Dental Association Council on Scientific Affairs

Douglas A. Young, DDS, EdD, MBA, MS; Brian B. Nový, DDS; Gregory G. Zeller, DDS, MS; Robert Hale, DDS; Thomas C. Hart, DDS, PhD; Edmond L. Truelove, DDS, MSD; American Dental Association Council on Scientific Affairs

Dental caries remains a common chronic disease and, in the absence of treatment, it may progress until the tooth is destroyed. Despite advances in restorative materials and the implementation of various preventive approaches, more than 90% of adults in the United States have experienced dental caries before 30 years of age.^{1,2}

Dental caries is a multifactorial disease involving many complex risk and protective factors.³ The clinical presentation of caries disease is a caries lesion; the severity of the disease and of individual caries lesions is the result of complex personal, biological, behavioral, and environmental factors. Some factors are protective, such as the presence of fluoride in the biofilm, whereas others lead to hard tissue destruction, such as lower plaque pH.⁴⁻⁶ Caries risk assessment is the organized process of evaluating these protective and pathogenic factors and provides the foundation⁷⁻⁹ for selecting treatment interventions.

The dental profession continues to implement a more interceptive nonsurgical therapeutic model to prevent, treat, and reverse caries lesions, particularly in the early stages. Despite progress, the profession still



This article has an accompanying online continuing education activity available at: <http://jada.ada.org/ce/home>. Copyright © 2015 American Dental Association. All rights reserved.

ABSTRACT

Background. The caries lesion, the most commonly observed sign of dental caries disease, is the cumulative result of an imbalance in the dynamic demineralization and remineralization process that causes a net mineral loss over time. A classification system to categorize the location, site of origin, extent, and when possible, activity level of caries lesions consistently over time is necessary to determine which clinical treatments and therapeutic interventions are appropriate to control and treat these lesions.

Methods. In 2008, the American Dental Association (ADA) convened a group of experts to develop an easy-to-implement caries classification system. The ADA Council on Scientific Affairs subsequently compiled information from these discussions to create the ADA Caries Classification System (CCS) presented in this article.

Conclusions. The ADA CCS offers clinicians the capability to capture the spectrum of caries disease presentations ranging from clinically unaffected (sound) tooth structure to noncavitated initial lesions to extensively cavitated advanced lesions. The ADA CCS supports a broad range of clinical management options necessary to treat both noncavitated and cavitated caries lesions.

Practical Implications. The ADA CCS is available for implementation in clinical practice to evaluate its usability, reliability, and validity. Feedback from clinical practitioners and researchers will allow system improvement. Use of the ADA CCS will offer standardized data that can be used to improve the scientific rationale for the treatment of all stages of caries disease.

Key Words. Caries classification system; caries lesion classification; caries location; caries extent; caries activity; caries management

JADA 2015;146(2):79-86


<http://dx.doi.org/10.1016/j.jada.2014.11.018>



The Importance of Cavitation

Bacteria are too big to penetrate the enamel





















When do you drill a tooth?



When it is fully cavitated

JADA
Feb 2015

TABLE 2

AMERICAN DENTAL ASSOCIATION CARIES CLASSIFICATION SYSTEM							
	Sound	Initial		Moderate	Advanced		
Clinical Presentation	No clinically detectable lesion. Dental hard tissue appears normal in color, translucency, and gloss.	Earliest clinically detectable lesion compatible with mild demineralization. Lesion limited to enamel or to shallow demineralization of cementum/dentin. Mildest forms are detectable only after drying. When established and active, lesions may be white or brown and enamel has lost its normal gloss.		Visible signs of enamel breakdown or signs the dentin is moderately demineralized.	Enamel is fully cavitated and dentin is exposed. Dentin lesion is deeply/severely demineralized.		
Other Labels	No surface change or adequately restored	Visually uncavitated		Established, early cavitated, shallow cavitation, microcavitation	Spread/disseminated, late cavitated, deep cavitation		
Infected Dentin	None	Unlikely		Possible	Present		
Appearance of Occlusal Surfaces (Pit and Fissure)*, †	ICDAS 0 	ICDAS 1 	ICDAS 2 	ICDAS 3 	ICDAS 4 	ICDAS 5 	ICDAS 6 
Accessible Smooth Surfaces, Including Cervical and Root‡							
Radiographic Presentation of the Approximal Surface§	 E0 [¶] or R0 [¶] No radiolucency	 E1 [¶] or RA1 [¶]	 E2 [¶] or RA2 [¶]	 D1 [¶] or RA3 [¶]	 D2 [¶] or RB4 [¶] Radiolucency extends into the middle one-third of the dentin	 D3 [¶] or RC5 [¶] Radiolucency extends into the inner one-third of the dentin	

Secondary

10

teeth

%

Primary

3

teeth

%

40

%

28

%

Use minimally invasive GI restoratives for moderate (D2) and advanced (D3) caries lesions only

2018 ADA Clinical Practice Guidelines



Cover Story

Evidence-based clinical practice guideline on nonrestorative treatments for carious lesions

A report from the American Dental Association

Rebecca L. Slayton, DDS, PhD; Olivia Urquhart, MPH; Marcelo W.B. Araujo, DDS, MS, PhD; Margherita Fontana, DDS, PhD; Sandra Guzmán-Armstrong, DDS, MS; Marcelle M. Nascimento, DDS, MS, PhD; Brian B. Nový, DDS; Norman Tinanoff, DDS, MS; Robert J. Weyant, DMD, DrPH; Mark S. Wolff, DDS, PhD; Douglas A. Young, DDS, EdD, MS, MBA; Domenick T. Zero, DDS, MS; Malavika P. Tampi, MPH; Lauren Pilcher, MSPH; Laura Banfield, MLIS, MHSc; Alonso Carrasco-Labra, DDS, MSc

ABSTRACT

Background. An expert panel convened by the American Dental Association Council on Scientific Affairs and the Center for Evidence-Based Dentistry conducted a systematic review and formulated evidence-based clinical recommendations for the arrest or reversal of noncavitated and cavitated dental caries using nonrestorative treatments in children and adults.

Types of Studies Reviewed. The authors conducted a systematic search of the literature in MEDLINE and Embase via Ovid, Cochrane CENTRAL, and Cochrane database of systematic reviews to identify randomized controlled trials reporting on nonrestorative treatments for noncavitated and cavitated carious lesions. The authors used the Grading of Recommendations Assessment, Development and Evaluation approach to assess the certainty in the evidence and move from the evidence to the decisions.

Results. The expert panel formulated 11 clinical recommendations, each specific to lesion type, tooth surface, and dentition. Of the most effective interventions, the panel provided recommendations for the use of 38% silver diamine fluoride, sealants, 5% sodium fluoride varnish, 1.23% acidulated phosphate fluoride gel, and 5,000 parts per million fluoride (1.1% sodium fluoride) toothpaste or gel, among others. The panel also provided a recommendation against the use of 10% casein phosphopeptide–amorphous calcium phosphate.

Conclusions and Practical Implications. Although the recommended interventions are often used for caries prevention, or in conjunction with restorative treatment options, these approaches have shown to be effective in arresting or reversing carious lesions. Clinicians are encouraged to prioritize use of these interventions based on effectiveness, safety, and feasibility.

SELECTIVE CARRIES REMOVAL

July 2023
JADA

[Clinical Practice Guideline](#)

[Cover Story](#)

Evidence-based clinical practice guideline on restorative treatments for caries lesions

A report from the American Dental Association



"More conservative carious tissue removal (CTR) approaches may decrease the risk of adverse events."



Types of Studies Reviewed. The authors searched for systematic reviews comparing carious tissue removal (CTR) approaches in Ovid MEDLINE, Embase, Cochrane Database of Systematic Reviews, and Trip Medical Database. The authors also conducted a systematic search for randomized controlled trials comparing direct restorative materials in Ovid MEDLINE, Embase, Cochrane Central Register of Controlled Trials, [ClinicalTrials.gov](#), and the World Health Organization International Clinical Trials Registry Platform. The authors used the Grading of Recommendations Assessment, Development, and Evaluation approach to assess the certainty of the evidence and formulate recommendations.

Results. The panel formulated 16 recommendations and good practice statements: 4 on CTR approaches specific to lesion depth and 12 on direct restorative materials specific to tooth location and surfaces involved. The panel conditionally recommended for the use of conservative CTR approaches, especially for advanced lesions. Although the panel conditionally recommended for the use of all direct restorative materials, they prioritized some materials over the use of others for certain clinical scenarios.



Practical Implications. The evidence suggests that more conservative CTR approaches may decrease the risk of adverse effects. All included direct restorative materials may be effective in treating moderate and advanced caries lesions on vital, nonendodontically treated primary and permanent teeth.

Key Words. Evidence-based dentistry; clinical practice guideline; direct restorative materials; caries; general dentistry; pediatric dentistry; American Dental Association.

Selective Caries Removal

Know when to stop!

- Clean the perimeter of the lesion using the 2mm rule to **avoid pulp exposure**
- Rationale: demineralization precedes bacterial penetration of dentinal tubules, so removing most of the **soft & wet** dentin will remove most of the infected tissue
- Glass ionomer will **seal out nutrient sources** from entering. Thus no acid or demineralization will occur



Infected = Wet & Soft

Affected = Dry & Soft

Selective caries removal will leave a radiolucency on the radiograph



Chemical exchange between glass-ionomer restorations and residual carious dentine in permanent molars: An in vivo study

Hien C. Ngo^{a,*}, Graham Mount^a, John Mc Intyre^a, J. Tuisuva^b, R.J. Von Doussa^a

^aDental School, University of Adelaide, Adelaide, Australia

^bDental School, Department of Medicine, Colonial War Memorial Hospital, Suva, Fiji

ARTICLE INFO

Article history:

Received 6 December 2004

Accepted 18 December 2004

Keywords:

Glass-ionomer

Dentine remineralization

Clinical trial

ABSTRACT

Objective: To evaluate the remineralization of carious dentine following the restoration of an extensive lesion in a permanent molar with a high strength glass-ionomer cement (GIC). **Materials and methods:** Thirteen first permanent molars, which were scheduled for extraction because of the presence of extensive caries lesions, were selected for this study. They were first restored, according to the ART technique, using encapsulated Fuji IX_{CP}, which contains a strontium glass rather than the traditional calcium glass. The cavities were prepared with a clean enamel margin and minimal removal of the carious dentine around the walls. After a period of 1–3 months they were harvested and subsequently sectioned and examined using an electron probe microanalysis (EPMA) and scanning electron microscopy (SEM).

Results: EPMA demonstrated that both fluoride and strontium ions had penetrated deep into the underlying demineralized dentine. The only possible source of these ions was the GIC restoration.

Conclusion: The pattern of penetration of the fluoride and strontium ions into the dentine was consistent with a remineralization process.

© 2006 Published by Elsevier Ltd.

1. Introduction

Since the time of Dr. G.V. Black the profession has been taught to completely remove softened and discoloured dentine to eliminate infected tissue and create a hard foundation to support a proposed restoration. The suggested routine has been to remove all demineralized dentine, using aggressive hand instrumentation or a round bur, until sound, normal dentine formed the entire pulpal floor. The objective was to ensure the elimination of all remaining microorganisms thus eliminating a possible recurrence of caries. However, Lager et al. showed that this is not always successful and some microorganisms may remain even after

all softened dentine has been removed and the cavity treated with sodium hypochlorite.¹ The main risk with this traditional approach is the possible accidental exposure of the pulp, particularly in young patients, where the rate of pulp exposure following excavation of large carious lesions in permanent molars has been rated at 40%.²

A step-wise excavation technique was introduced by Bodecker³ designed to decrease the risk of mechanical pulp exposure. Bodecker recommended partial removal of the soft demineralized dentine on the cavity floor followed by immediate restoration with a temporary material such as zinc oxide/eugenol. The transitional material was expected to remain for a brief period of weeks and then replaced with

Internal Remineralization

“Electron probe microanalysis demonstrated that both fluoride and strontium ions had penetrated deep into underlying demineralized dentin. The pattern was consistent with remineralization. The only source of these ions was the glass ionomer restoration.”

* Corresponding author. Tel.: +61 8 8303 5256; fax: +61 8 8303 3444.

E-mail address: hien.ngo@adelaide.edu.au (H.C. Ngo).


0300-5712/\$ – see front matter © 2006 Published by Elsevier Ltd.

doi:10.1016/j.jdent.2005.12.012


Studies by Dr. Hien Ngo

The "Chemically Fused Layer on Enamel"

Internal Remineralization is Acid Resistant



Chemically Fused Zone on Dentin



SEM technique for examining the glass-ionomer cement chemical fusion zone. Milicich G. Journal of Microscopy, Vol. 217, Pt 1 January 2005, pp. 44–48

Practical implementation of caries disease management in the era of CAMBRA: a guide for the general dentist

Douglas A. Young, DDS, EdD, MBA, MS ■ Ryan L. Quock, DDS ■ Susie Goolsby, DDS, MSHA
 Martha McComas, RDH, MS ■ Steven Powell, DDS ■ Gregory G. Zeller, DDS, MS
 Allen Wong, DDS, EdD ■ Brian B. Nový, DDS

Although the philosophy of caries management by risk assessment (CAMBRA) is widely accepted as an essential part of comprehensive disease management, its implementation in daily clinical practice has historically been hindered by a lack of concise guidance that bridges the gap between scholarly research and the logistics of patient care. With the aim of guiding comprehensive disease management for dental caries by integrating available evidence-based systems such as the American Dental Association Caries Classification System (ADA CCS) and caries risk assessment (CRA), this review presents clinical strategies in a functional format. Clinical caries management strategies that integrate the ADA CCS with CRA were compiled and evaluated by clinicians representing a variety of practice settings. Examples of clinical protocols were selected to illustrate the use of CRA at the patient level and ADA CCS lesion classification at the tooth level to comprehensively manage caries lesions in a high-carries-risk population. Comprehensive caries management combines the clinical use of CRA and caries lesion classification to guide clinical decisions and to support tooth preservation in high-risk populations. These strategies should be individualized and fulfill evidence-based clinical practice criteria with oral health as the primary outcome. Person-centered comprehensive caries management integrating CRA and classification of caries lesions has the potential to improve patient outcomes by addressing the underlying causes of caries disease. Documenting and measuring the processes at the tooth, patient, and practice levels are essential for generating outcomes that will drive future approaches aimed at improving patient health.

Received: February 16, 2025
 Revised: May 25, 2025
 Accepted: June 30, 2025

Keywords: caries classification system, caries management, caries risk assessment, caries risk management, clinical protocols, dental caries, oral health

Published with permission of the Academy of General Dentistry. © Copyright 2025 by the Academy of General Dentistry. All rights reserved. For printed and electronic reprints of this article for distribution, please contact scsreprints@sheridan.com.



















The respect for comprehensive disease management as a whole, and caries management by risk assessment (CAMBRA) as a leading philosophy, has continued to increase among payers, educators, clinicians, organized oral healthcare, and the US Health Resources and Services Administration. Yet in an era of evidence-based clinical decision-making, the adoption of these concepts by more clinicians seems to be stalled by a lack of practical translational strategies that are urgently needed for the improvement of patient health. This is not to say that the underlying science is lacking or that research has failed to show any positive health outcomes from these approaches. Instead, the situation highlights the hurdles that make it difficult to evaluate and adopt the vast array of available information in the form of effective clinical protocols. Given that dental caries is one of the most common diseases treated in dentistry, it could be argued that every general dentist functions as a practical cariologist, requiring actionable information. This summary of clinical cariology based on the philosophy of CAMBRA is meant to provide guidance for general dentists to implement disease management concepts in daily practice and improve the oral health of the patients they treat.

Terminology

The term *caries* by itself is open to several interpretations, and *caries* is colloquially used to describe both the disease and the resulting symptomatic lesion; this may create confusing variations in diagnosis and treatment. Adoption of high-level consensus terminology, such as that described jointly by the International Association for Dental Research Cariology Research Group, the European Organization for Caries Research, and the International Caries Consensus Collaboration is recommended for more precise diagnosis, treatment, and communication.¹⁻³ For example, the term *dental caries* is used to describe the actual disease process, while *caries lesion* describes the clinical manifestation of the disease; the 2 terms are closely related yet distinct.


Caries management is defined as "...actions taken to interfere with mineral loss at all stages of the caries disease, including nonoperative and operative interventions/treatments."¹ The term *nonoperative* refers to treating the lesion without tooth structure removal by employing the following strategies: surface protection (including barriers or sealants), biofilm modification, remineralization via chemotherapeutics, saliva supplementation, and behavioral modification, including changes in diet and lifestyle.

Table 5. Example hard tissue preservation strategy based on the ADA CCS category using current tooth-preserving (tooth-level) and behavior management (patient-level) strategies.

Characteristics	Sound	Initial	Moderate	Advanced			
Infected dentin (full cavitation)	None	Unlikely	Possible	Present			
Occlusal							
Clinical presentation	ICDAS 0 	ICDAS 1 	ICDAS 2 	ICDAS 3 	ICDAS 4 	ICDAS 5 	ICDAS 6 
Insidious caries lesion; few lesions, slow progression	Consider GIC sealant for primary prevention (no preparation)	Place GIC sealant (no preparation)	Place GIC sealant (no preparation)	Place MID GIC restoration (selective caries removal)			
Rampant or aggressive caries lesion; multiple lesions, aggressive progression	Place GIC sealant (no preparation)	Arrest with SDF; place GIC sealant (no preparation)	Arrest with SDF; place GIC sealant (no preparation)	Consider SDF; place MID GIC restoration (selective caries removal)			
Facial/lingual							
Clinical presentation	ICDAS 0 	ICDAS 1 	ICDAS 2 	ICDAS 3 	ICDAS 4 	ICDAS 5 	ICDAS 6 
Insidious caries lesion; few lesions, slow progression	Follow CRM protocol for primary prevention	Arrest and remineralize using risk-based CRM; consider SDF	Arrest and remineralize using risk-based CRM; consider SDF surface protection (no preparation) for deeper microcavitations and/or SDF per patient wants	Place MID GIC restoration (selective caries removal)			
Rampant or aggressive caries lesion; multiple lesions, aggressive progression	Follow CRM protocol for primary prevention	Arrest and remineralize using risk-based CRM; consider SDF	Arrest and remineralize using risk-based CRM; consider GIC surface protection (no preparation) and/or SDF	Consider SDF; place GIC restoration (selective caries removal)			
Approximal							
Clinical presentation	D1* or R0* 	E1* or RA1* E2* or RA2* D1* or RA1* 	D2* or RB4* 	D3* or RC5* 			
Difficult to assess insidious vs aggressive progression	Follow CRM protocol for primary prevention	Arrest and remineralize using risk-based CRM; consider SDF	Arrest and remineralize using risk-based CRM; consider SDF orthodontic separation, surface protection with no preparation, and/or MID restoration with selective caries removal with or without SDF	Place MID GIC restoration (selective caries removal) with or without SDF			
CRM protocol: topical fluoride; pH neutralization; biofilm modification; saliva supplementation with calcium and phosphate based on caries risk; and behavioral modification (healthy diet and lifestyle choices and good oral hygiene).							
Abbreviations: ADA CCS, American Dental Association Caries Classification System; CRM, caries risk management; GIC, glass ionomer cement; ICDAS, International Caries Detection and Assessment System; MID, minimally invasive dentistry; SDF, silver diamine fluoride.							
Images reproduced from Young DA, Nový BB, Zeller GG, et al. The American Dental Association Caries Classification System for Clinical Practice: a report of the American Dental Association Council on Scientific Affairs. <i>J Am Dent Assoc.</i> 2015;146(2):79-86. doi:10.1016/j.adaj.2014.11.018. Copyright ©2015 American Dental Association. All rights reserved. Reprinted by permission.							

Agd.org/generaldentistry; Young DA, Quock RL, Goolsby S, et al. Practical implementation of caries disease management in the era of CAMBRA: a guide for the general dentist. *Gen Dent.* 2025;73(6):56-66.

SMART Textbook



PDF now available for free in English and Spanish from the following website

www.smartoralhealth.com